



Product data sheet REF 23.40.20

25.08.2026

General data			
article number	23.40.20	Kocher retr 20x38mm 23cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	45918	GTIN	04046826031077
Risk category	1r (EU MDR)	UMDNS-Code	15-242
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;03;aj;5L
Net weight (kg)	0,100	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK12-1



Medicon media	450.00.03 (general catalogue the art of surgery)
---------------	--

Complete picture

Detail images



size	20x38mm 23cm	Coating	
------	--------------	---------	--

Storage conditions

No special requirements

Materials

	Material number (DIN designation)
Main part	1.4021 (X20Cr13)