



Product data sheet REF 23.40.55

25.08.2026

General data			
article number	23.40.55	Kocher retr 55x80mm 25.5cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	45918	GTIN	04046826031107
Risk category	1r (EU MDR)	UMDNS-Code	15-242
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;03;aj;5L
Net weight (kg)	0,193	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK24-1
			
Medicon media	450.00.03 (general catalogue the art of surgery)		
Complete picture		Detail images	
			
size	55x80mm 25.5cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		