



Product data sheet REF 23.54.31

11.12.2025

| General data | | | |
|-----------------|-------------|----------------------------------|-------------------------|
| article number | 23.54.31 | Doyen retr/speculum 35x55mm 24cm | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 45918 | GTIN | 04046826031435 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 13-666 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R004;03;aj;5L |
| Net weight (kg) | 0,180 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | MPK20-1 |



| | |
|---------------|--|
| Medicon media | 450.00.03 (general catalogue the art of surgery) |
|---------------|--|

Complete picture

Detail images



| | | | |
|------|--------------|---------|--|
| size | 35x55mm 24cm | Coating | |
|------|--------------|---------|--|

Storage conditions

No special requirements

Materials

| | Material number (DIN designation) |
|-----------|-----------------------------------|
| Main part | 1.4301 (X5CrNi18-10) |
| Handle | 1.4301 (X5CrNi18-10) |