



Product data sheet REF 23.54.34

03.06.2026

General data			
article number	23.54.34	Doyen retr/speculum 35x150mm 24cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	45918	GTIN	04046826031466
Risk category	1r (EU MDR)	UMDNS-Code	13-666
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;03;aj;5L
Net weight (kg)	0,255	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK24-1
			
Medicon media	450.00.03 (general catalogue the art of surgery)		
Complete picture		Detail images	
			
size	150x35mm 24cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4301 (X5CrNi18-10)		
Handle	1.4301 (X5CrNi18-10)		