



## Product data sheet REF 23.54.64

03.06.2026

| General data                                                                                                                                                                  |                                                  |                                   |                         |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------|-----------------------------------|-------------------------|
| article number                                                                                                                                                                | 23.54.64                                         | Doyen retr/speculum 60x155mm 24cm |                         |
| Unit                                                                                                                                                                          | piece                                            | Customs code (HS code)            | 90189084                |
| GMDN                                                                                                                                                                          | 45918                                            | GTIN                              | 04046826031541          |
| Risk category                                                                                                                                                                 | 1r (EU MDR)                                      | UMDNS-Code                        | 13-666                  |
| Classification                                                                                                                                                                | 6 (EU MDR)                                       | Base UDI-DI                       | 4046826;B;R004;03;aj;5L |
| Net weight (kg)                                                                                                                                                               | 0,290                                            | Operating instructions            | GATDR00X                |
| Latex free                                                                                                                                                                    | yes                                              | Preparation instructions          | MPK24-1                 |
| <div></div> |                                                  |                                   |                         |
| Medicon media                                                                                                                                                                 | 450.00.03 (general catalogue the art of surgery) |                                   |                         |
| Complete picture                                                                                                                                                              |                                                  | Detail images                     |                         |
|                                                                                             |                                                  |                                   |                         |
| size                                                                                                                                                                          | 60x155mm 24cm                                    | Coating                           |                         |
| Storage conditions                                                                                                                                                            |                                                  |                                   |                         |
| No special requirements                                                                                                                                                       |                                                  |                                   |                         |
| Materials                                                                                                                                                                     |                                                  |                                   |                         |
|                                                                                                                                                                               | Material number (DIN designation)                |                                   |                         |
| Main part<br>Handle                                                                                                                                                           | 1.4021 (X20Cr13)                                 |                                   |                         |
|                                                                                                                                                                               | 1.4301 (X5CrNi18-10)                             |                                   |                         |