



Product data sheet REF 23.72.17

03.06.2026

General data			
article number	23.72.17	Deaver Baby retr 17mm 20cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	45918	GTIN	04046826031978
Risk category	1r (EU MDR)	UMDNS-Code	13-649
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;03;aj;5L
Net weight (kg)	0,050	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK12-1
			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.53 (catalogue cardio)		
Complete picture		Detail images	
size	17mm 20cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		