



Product data sheet REF 23.73.50

03.06.2026

General data			
article number	23.73.50	Deaver retr 50mm 32cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	45918	GTIN	04046826032012
Risk category	1r (EU MDR)	UMDNS-Code	13-649
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;03;aj;5L
Net weight (kg)	0,259	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK24-1
			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.53 (catalogue cardio)		
Complete picture		Detail images	
			
size	50mm 32cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part Handle	1.4021 (X20Cr13) 1.4301 (X5CrNi18-10)		