



Product data sheet REF 23.75.06

03.06.2026

General data			
article number	23.75.06	Deaver retr 25mm 35cm fig 2.5	
Unit	piece	Customs code (HS code)	90189084
GMDN	45918	GTIN	04046826032081
Risk category	1r (EU MDR)	UMDNS-Code	13-649
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;03;aj;5L
Net weight (kg)	0,131	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK20-1



Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.53 (catalogue cardio)
---------------	---

Complete picture

Detail images



size	25mm 35cm	Coating	
------	-----------	---------	--

Storage conditions

No special requirements

Materials

	Material number (DIN designation)
Main part	1.4021 (X20Cr13)