



Product data sheet REF 23.75.09

03.06.2026

General data			
article number	23.75.09	Deaver retr 75mm 30cm fig 5	
Unit	piece	Customs code (HS code)	90189084
GMDN	45918	GTIN	04046826032111
Risk category	1r (EU MDR)	UMDNS-Code	13-649
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;03;aj;5L
Net weight (kg)	0,233	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK24-1



Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.53 (catalogue cardio)
---------------	---

Complete picture

Detail images



size	75mm 30cm	Coating	
------	-----------	---------	--

Storage conditions

No special requirements

Materials

	Material number (DIN designation)
Main part	1.4021 (X20Cr13)