



## Product data sheet REF 24.10.70

25.08.2026

| General data                                                                                  |                                                  |                                      |                    |
|-----------------------------------------------------------------------------------------------|--------------------------------------------------|--------------------------------------|--------------------|
| article number                                                                                | 24.10.70                                         | Logan (Goldstein) retr 3x3t sh 3.5cm |                    |
| Unit                                                                                          | piece                                            | Customs code (HS code)               | 90189084           |
| GMDN                                                                                          | 45182                                            | GTIN                                 | 04046826032517     |
| Risk category                                                                                 | 2a (EU MDD)                                      | UMDNS-Code                           | 13-390             |
| Classification                                                                                | 6 (EU MDD)                                       | Base UDI-DI                          | 4046826;A002;03;P6 |
| Net weight (kg)                                                                               | 0,007                                            | Operating instructions               |                    |
| Latex free                                                                                    | yes                                              | Preparation instructions             | MPK02              |
| <div></div>  |                                                  |                                      |                    |
| Medicon media                                                                                 | 450.00.03 (general catalogue the art of surgery) |                                      |                    |
| Complete picture                                                                              |                                                  | Detail images                        |                    |
| <div></div> |                                                  |                                      |                    |
| size                                                                                          | 3.5cm                                            | Coating                              |                    |
| Storage conditions                                                                            |                                                  |                                      |                    |
| No special requirements                                                                       |                                                  |                                      |                    |
| Materials                                                                                     |                                                  |                                      |                    |
|                                                                                               | Material number (DIN designation)                |                                      |                    |
| Main part                                                                                     | 1.4021 (X20Cr13)                                 |                                      |                    |
| Others                                                                                        | 1.4104 (X14CrMoS17)                              |                                      |                    |