



Product data sheet REF 24.23.01

11.12.2025

General data			
article number	24.23.01	Funk tissue retractor 90mm left	
Unit	piece	Customs code (HS code)	90189084
GMDN	45918	GTIN	04046826457112
Risk category	1r (EU MDR)	UMDNS-Code	15-218
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;02;aj;59
Net weight (kg)	0,000	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK23-1
 			
Medicon media	453.08.00 (Presto Brochure)		
Complete picture		Detail images	
			
size	90mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		
Screw	1.4305 (X8CrNiS18-9)		
Others	1.4310 (X10CrNi18-8)		