



Product data sheet REF 24.26.06

25.08.2026

General data			
article number	24.26.06	Wigand retr 10.5cm rt 3 pins lf 3t	
Unit	piece	Customs code (HS code)	90189084
GMDN	45182	GTIN	04046826409791
Risk category	2a (EU MDD)	UMDNS-Code	13-390
Classification	6 (EU MDD)	Base UDI-DI	4046826;A002;03;P6
Net weight (kg)	0,044	Operating instructions	
Latex free	yes	Preparation instructions	MPK03
			
Medicon media	450.00.37 (catalogue ENT)		
Complete picture		Detail images	
			
size	10.5cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		