



Product data sheet REF 24.34.45

11.12.2025

General data			
article number	24.34.45	Beckmann Adson retr 44mm sh 31cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	45182	GTIN	04046826033125
Risk category	2a (EU MDD)	UMDNS-Code	13-390
Classification	7 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,259	Operating instructions	
Latex free	yes	Preparation instructions	MPK03
			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.53 (catalogue cardio)		
Complete picture		Detail images	
size	44mm 31cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		
Screw	1.4305 (X8CrNiS18-9)		
Others	1.4310 (X10CrNi18-8)		