



## Product data sheet REF 24.70.30

11.12.2025

| General data                                                                                  |                                                  |                                    |                |
|-----------------------------------------------------------------------------------------------|--------------------------------------------------|------------------------------------|----------------|
| article number                                                                                | 24.70.30                                         | Balfour Medicon abdom retr w ratch |                |
| Unit                                                                                          | piece                                            | Customs code (HS code)             | 90189084       |
| GMDN                                                                                          | 45182                                            | GTIN                               | 04046826033903 |
| Risk category                                                                                 | 2a (EU MDD)                                      | UMDNS-Code                         | 13-390         |
| Classification                                                                                | 7 (EU MDD)                                       | Base UDI-DI                        |                |
| Net weight (kg)                                                                               | 1,150                                            | Operating instructions             |                |
| Latex free                                                                                    | yes                                              | Preparation instructions           | MPK02          |
| <div></div>  |                                                  |                                    |                |
| Medicon media                                                                                 | 450.00.03 (general catalogue the art of surgery) |                                    |                |
| Complete picture                                                                              |                                                  | Detail images                      |                |
| <div></div> |                                                  |                                    |                |
|                                                                                               |                                                  |                                    |                |
| Storage conditions                                                                            |                                                  |                                    |                |
| No special requirements                                                                       |                                                  |                                    |                |
| Materials                                                                                     |                                                  |                                    |                |
|                                                                                               | Material number (DIN designation)                |                                    |                |
| Main part                                                                                     | 1.4301 (X5CrNi18-10)                             |                                    |                |
| Screw                                                                                         | 1.4305 (X8CrNiS18-9)                             |                                    |                |