





Product data sheet REF 24.70.33

11.12.2025

| General data | | | |
|--|--|--------------------------|----------------|
| article number | 24.70.33 | lateral blade 90mm pair | |
| Unit | box/2 | Customs code (HS code) | 90189084 |
| GMDN | 45182 | GTIN | 04046826033934 |
| Risk category | 2a (EU MDD) | UMDNS-Code | 13-390 |
| Classification | 7 (EU MDD) | Base UDI-DI | |
| Net weight (kg) | 0,065 | Operating instructions | |
| Latex free | yes | Preparation instructions | MPK02 |
|  | | | |
| Medicon media | 450.00.03 (general catalogue the art of surgery) | | |
| Complete picture | | Detail images | |
|  | | | |
| size | 90mm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part Screw | 1.4021 (X20Cr13) 1.4305 (X8CrNiS18-9) | | |