








## Product data sheet REF 26.24.05

11.12.2025

General data			
article number	26.24.05	Williger elevator 5.5mm 17cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	35048	GTIN	04046826034948
Risk category	1r (EU MDR)	UMDNS-Code	15-218
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;05;aj;6A
Net weight (kg)	0,029	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK03-1
 			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.37 (catalogue ENT ) 450.00.40 (Katalog Nr 5 Dental ) 450.00.63 (catalogue Orthopaedics )		
Complete picture		Detail images	
  			
size	5.5mm 17cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4024 (X15Cr13)		