







Product data sheet REF 26.27.70

11.12.2025

General data			
article number	26.27.70	Hoen elevator 20mm 19cm semi sh	
Unit	piece	Customs code (HS code)	90189084
GMDN	35048	GTIN	04046826035020
Risk category	1r (EU MDR)	UMDNS-Code	15-218
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;05;aj;6A
Net weight (kg)	0,096	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK12-1
 			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.63 (catalogue Orthopaedics)		
Complete picture		Detail images	
 			
size	20mm 19cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		