







## Product data sheet REF 26.73.09

11.12.2025

| General data   |   |                          |                         |
|--|---|--------------------------|-------------------------|
| article number   | 26.73.09  | raspatory 9mm 28cm       |                         |
| Unit   | piece   | Customs code (HS code)   | 90189084                |
| GMDN   | 35559   | GTIN                     | 04046826036119          |
| Risk category  | 1r (EU MDR)   | UMDNS-Code               | 11-084                  |
| Classification   | 6 (EU MDR)  | Base UDI-DI              | 4046826;B;R003;03;ai;4W |
| Net weight (kg)  | 0,150   | Operating instructions   | GATDR00X                |
| Latex free   | yes   | Preparation instructions | MPK20-1                 |
|      |   |                          |                         |
| Medicon media  | 450.00.03 (general catalogue the art of surgery)<br>450.00.63 (catalogue Orthopaedics ) |                          |                         |
| Complete picture   |   | Detail images            |                         |
|   |   |                          |                         |
| size   | 9mm 28cm  | Coating                  |                         |
| Storage conditions   |   |                          |                         |
| No special requirements  |   |                          |                         |
| Materials  |   |                          |                         |
|  | Material number (DIN designation)   |                          |                         |
| Main part<br>Handle  | 1.4021 (X20Cr13)<br>1.4301 (X5CrNi18-10)  |                          |                         |