





## Product data sheet REF 27.10.73

11.12.2025

| General data  |   |                          |                |
|---|---|--------------------------|----------------|
| article number  | 27.10.73  | wrench f 27.10.70        |                |
| Unit  | piece   | Customs code (HS code)   | 90189084       |
| GMDN  | 46976   | GTIN                     | 04046826036560 |
| Risk category   | 1 (EU MDD)  | UMDNS-Code               | 15-229         |
| Classification  | 1 (EU MDD)  | Base UDI-DI              |                |
| Net weight (kg)   | 0,020   | Operating instructions   |                |
| Latex free  | yes   | Preparation instructions | MPK01          |
|    |   |                          |                |
| Medicon media   | 450.00.03 (general catalogue the art of surgery)<br>450.00.63 (catalogue Orthopaedics ) |                          |                |
| Complete picture  |   | Detail images            |                |
|  |   |                          |                |
| Storage conditions  |   |                          |                |
| No special requirements   |   |                          |                |
| Materials   |   |                          |                |
|   | Material number (DIN designation)   |                          |                |
| Main part   | 1.4021 (X20Cr13)  |                          |                |