



Product data sheet REF 277.53.90

11.12.2025

| General data | | | |
|--|---|--|--|
| article number | 277.53.90 | Bl-Line SC Converse Daniel scrs ang 9.5cm sh | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 13495 | GTIN | 04046826445058 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 13-500 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R001;01;ag;2U |
| Net weight (kg) | 0,000 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | MPK23-1 |
|   | | | |
| Medicon media | 450.00.60 (catalogue plastic surgery) 453.30.00 (daniel - Advanced Rhinoplasty Instrument Set Brochure) | | |
| Complete picture | | Detail images | |
|   | | | |
| size | 9.5cm | Coating | Artificially created PVD layer of titanium aluminum carbon nitrite e.g. Magic Cut scissors |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part Screw | 1.4117 (X38CrMoV15) 1.4305 (X8CrNiS18-9) | | |