



## Product data sheet REF 28.42.08

11.12.2025

| General data    |             |                          |                         |
|-----------------|-------------|--------------------------|-------------------------|
| article number  | 28.42.08    | Stille gouge 8mm 16cm    |                         |
| Unit            | piece       | Customs code (HS code)   | 90189084                |
| GMDN            | 10825       | GTIN                     | 04046826038458          |
| Risk category   | 1r (EU MDR) | UMDNS-Code               | 10-825                  |
| Classification  | 6 (EU MDR)  | Base UDI-DI              | 4046826;B;R003;01;ai;48 |
| Net weight (kg) | 0,066       | Operating instructions   | GATDR00X                |
| Latex free      | yes         | Preparation instructions | MPK12-1                 |



Medicon media 450.00.03 (general catalogue the art of surgery)  
450.00.63 (catalogue Orthopaedics )

## Complete picture

## Detail images



|      |          |         |  |
|------|----------|---------|--|
| size | 8mm 16cm | Coating |  |
|------|----------|---------|--|

## Storage conditions

No special requirements

## Materials

|           | Material number (DIN designation) |
|-----------|-----------------------------------|
| Main part | 1.4034 (X46Cr13)                  |