



## Product data sheet REF 28.65.20P

11.12.2025

General data			
article number	28.65.20P	Lexer gouge 20mm 22cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	10825	GTIN	04046826432171
Risk category	1r (EU MDR)	UMDNS-Code	10-825
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;01;ai;48
Net weight (kg)	0,000	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK20-1



Medicon media	450.00.63 (catalogue Orthopaedics )
---------------	-------------------------------------

## Complete picture

## Detail images



size	20mm 22cm	Coating	
------	-----------	---------	--

## Storage conditions

No special requirements

## Materials

	Material number (DIN designation)
Main part Handle	1.4034 (X46Cr13) PROPYLUX (Polypropylen Homopolymere)