







Product data sheet REF 28.67.04

11.12.2025

| General data | | | |
|--|---|--------------------------|-------------------------|
| article number | 28.67.04 | Dahmen chisel 4mm 30cm | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 10825 | GTIN | 04046826039912 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 10-825 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R003;01;ai;48 |
| Net weight (kg) | 0,140 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | |
|   | | | |
| Medicon media | 450.00.03 (general catalogue the art of surgery) 450.00.63 (catalogue Orthopaedics) | | |
| Complete picture | | Detail images | |
|   | | | |
| size | 4mm 30cm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part Handle | 1.4034 (X46Cr13) HGW_2082 (PF_CC_201) | | |