



Product data sheet REF 28.67.52P

11.12.2025

General data			
article number	28.67.52P	Dahmen chisel 22mm 30cm ang	
Unit	piece	Customs code (HS code)	90189084
GMDN	10825	GTIN	04046826432690
Risk category	1r (EU MDR)	UMDNS-Code	10-825
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;01;ai;48
Net weight (kg)	0,000	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK20-1
article number	28.67.52P	Dahmen chisel 22mm 30cm ang	
Unit	piece	Customs code (HS code)	90189084
GMDN	10825	GTIN	04046826432690
Risk category	1r (EU MDR)	UMDNS-Code	10-825
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;01;ai;48
Net weight (kg)	0,000	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK20-1




Medicon media	450.00.63 (catalogue Orthopaedics)
---------------	-------------------------------------

Medicon media	450.00.63 (catalogue Orthopaedics)
---------------	-------------------------------------

Complete picture

Detail images





size	22mm 30cm	Coating	
size	22mm 30cm	Coating	

Storage conditions

No special requirements

No special requirements

Materials

	Material number (DIN designation)
Main part Handle	1.4034 (X46Cr13) PROPYLUX (Polypropylen Homopolymere)
Main part Handle	1.4034 (X46Cr13) PROPYLUX (Polypropylen Homopolymere)