



Product data sheet REF 28.67.64P

11.12.2025

| General data | | | |
|-----------------|-------------|----------------------------|-------------------------|
| article number | 28.67.64P | Dahmen chisel 4mm 30cm ang | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 10825 | GTIN | 04046826431730 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 10-825 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R003;01;ai;48 |
| Net weight (kg) | 0,000 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | MPK20-1 |



| | |
|---------------|-------------------------------------|
| Medicon media | 450.00.63 (catalogue Orthopaedics) |
|---------------|-------------------------------------|

Complete picture

Detail images



| | | | |
|------|----------|---------|--|
| size | 4mm 30cm | Coating | |
|------|----------|---------|--|

Storage conditions

No special requirements

Materials

| | Material number (DIN designation) |
|---------------------|--|
| Main part Handle | 1.4034 (X46Cr13) PROPYLUX (Polypropylen Homopolymere) |