








Product data sheet REF 28.67.85P

11.12.2025

General data			
article number	28.67.85P	Dahmen chisel 25mm 30cm ang	
Unit	piece	Customs code (HS code)	90189084
GMDN	10825	GTIN	04046826431822
Risk category	1r (EU MDR)	UMDNS-Code	10-825
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;01;ai;48
Net weight (kg)	0,000	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK24-1
 			
Medicon media			
Complete picture		Detail images	
  			
size	25mm 30cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part Handle	1.4034 (X46Cr13) PROPYLUX (Polypropylen Homopolymere)		