







## Product data sheet REF 29.18.08

11.12.2025

| General data   |   |                                  |                         |
|--|---|----------------------------------|-------------------------|
| article number   | 29.18.08  | Uffenorde sh curette 45° rt 17cm |                         |
| Unit   | piece   | Customs code (HS code)           | 90189084                |
| GMDN   | 31335   | GTIN                             | 04046826040949          |
| Risk category  | 1r (EU MDR)   | UMDNS-Code                       | 11-084                  |
| Classification   | 6 (EU MDR)  | Base UDI-DI                      | 4046826;B;R003;02;ai;4K |
| Net weight (kg)  | 0,048   | Operating instructions           | GATDR00X                |
| Latex free   | yes   | Preparation instructions         | MPK12-1                 |
|      |   |                                  |                         |
| Medicon media  | 450.00.03 (general catalogue the art of surgery)<br>450.00.37 (catalogue ENT )<br>450.00.46 (catalogue sinuscopy )<br>450.00.63 (catalogue Orthopaedics ) |                                  |                         |
| Complete picture   |   | Detail images                    |                         |
|   |   |                                  |                         |
| size   | 45° 17cm  | Coating                          |                         |
| Storage conditions   |   |                                  |                         |
| No special requirements  |   |                                  |                         |
| Materials  |   |                                  |                         |
|  | Material number (DIN designation)   |                                  |                         |
| Main part<br>Handle  | 1.4021 (X20Cr13)<br>1.4301 (X5CrNi18-10)  |                                  |                         |