



Product data sheet REF 29.56.57

25.08.2026

General data			
article number	29.56.57	sh serrated curette 8x150mm 8x150mm	
Unit	piece	Customs code (HS code)	90189084
GMDN	31335	GTIN	04046826540463
Risk category	1r (EU MDR)	UMDNS-Code	11-084
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;02;ai;4K
Net weight (kg)	0,118	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK20-1
			
Medicon media			
Complete picture		Detail images	
			
size	8x150mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part Handle	1.4034 (X46Cr13) SILIKON (Poly(organo)siloxane)		