






Product data sheet REF 29.57.53P

11.12.2025

| General data | | | |
|---|--|-----------------------------|-------------------------|
| article number | 29.57.53P | Bruns sh curette fig 3 33cm | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 31335 | GTIN | 04046826432874 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 11-084 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R003;02;ai;4K |
| Net weight (kg) | 0,000 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | MPK20-1 |
|   | | | |
| Medicon media | | | |
| Complete picture | | Detail images | |
|  | | | |
| size | 33cm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part Handle | 1.4021 (X20Cr13) PROPYLUX (Polypropylen Homopolymere) | | |