








Product data sheet REF 29.57.70P

11.12.2025

General data			
article number	29.57.70P	curette serr 5mm 39cm strongly cvd	
Unit	piece	Customs code (HS code)	90189084
GMDN	31335	GTIN	04046826459550
Risk category	1r (EU MDR)	UMDNS-Code	13-508
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;02;ai;4K
Net weight (kg)	0,151	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	
 			
Medicon media			
Complete picture		Detail images	
  			
size	5mm 39cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part Handle	1.4021 (X20Cr13) PROPYLUX (Polypropylen Homopolymere)		