








## Product data sheet REF 29.57.78

11.12.2025

| General data  |   |                           |                         |
|---|---|---------------------------|-------------------------|
| article number  | 29.57.78  | curette serr 4mm 39cm cvd |                         |
| Unit  | piece   | Customs code (HS code)    | 90189084                |
| GMDN  | 31335   | GTIN                      | 04046826042042          |
| Risk category   | 1r (EU MDR)   | UMDNS-Code                | 13-508                  |
| Classification  | 6 (EU MDR)  | Base UDI-DI               | 4046826;B;R003;02;ai;4K |
| Net weight (kg)   | 0,200   | Operating instructions    | GATDR00X                |
| Latex free  | yes   | Preparation instructions  |                         |
|     |   |                           |                         |
| Medicon media   | 451.04.65 (brochure Thora MIC Thorakosk Mini Inv Surgery) |                           |                         |
| Complete picture  |   | Detail images             |                         |
|    |   |                           |                         |
| size  | 4mm 39cm  | Coating                   |                         |
| Storage conditions  |   |                           |                         |
| No special requirements   |   |                           |                         |
| Materials   |   |                           |                         |
|   | Material number (DIN designation)                         |                           |                         |
| Main part<br>Handle   | 1.4021 (X20Cr13)<br>HGW_2082 (PF_CC_201)                  |                           |                         |