







Product data sheet REF 30.07.15

11.12.2025

| General data | | | |
|---|---|-----------------------------|-------------------------|
| article number | 30.07.15 | bone rong 14cm strongly cvd | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 32853 | GTIN | 04046826042554 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 15-243 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R001;04;ag;3V |
| Net weight (kg) | 0,098 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | MPK23-1 |
|   | | | |
| Medicon media | 450.00.03 (general catalogue the art of surgery) 450.00.37 (catalogue ENT) 450.00.40 (Katalog Nr 5 Dental) 450.00.62 (catalogue Medicon instruments Oral and Maxillofacial Surgery) 450.00.63 (catalogue Orthopaedics) | | |
| Complete picture | | Detail images | |
|   | | | |
| size | 14cm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4021 (X20Cr13) | | |
| Others | 1.4021 (X20Cr13) | | |