







## Product data sheet REF 30.15.17

11.12.2025

| General data  |   |   |                         |
|---|---|---|-------------------------|
| article number  | 30.15.17  | Mead bone rong fig 4A 17cm  |                         |
| Unit  | piece   | Customs code (HS code)  | 90189084                |
| GMDN  | 32853   | GTIN  | 04046826042578          |
| Risk category   | 1r (EU MDR)   | UMDNS-Code  | 15-243                  |
| Classification  | 6 (EU MDR)  | Base UDI-DI   | 4046826;B;R001;04;ag;3V |
| Net weight (kg)   | 0,146   | Operating instructions  | GATDR00X                |
| Latex free  | yes   | Preparation instructions  | MPK23-1                 |
|   |   |   |                         |
| Medicon media   | 450.00.03 (general catalogue the art of surgery)<br>450.00.40 (Katalog Nr 5 Dental )<br>450.00.62 (catalogue Medicon instruments Oral and Maxillofacial Surgery)<br>450.00.63 (catalogue Orthopaedics ) |   |                         |
| Complete picture  |   | Detail images   |                         |
|   |   |  |                         |
| size  | 17cm  | Coating   |                         |
| Storage conditions  |   |   |                         |
| No special requirements   |   |   |                         |
| Materials   |   |   |                         |
|   | Material number (DIN designation)   |   |                         |
| Main part   | 1.4021 (X20Cr13)  |   |                         |