



Product data sheet REF 30.52.18

11.12.2025

| General data | | | |
|-----------------|-------------|--------------------------|-------------------------|
| article number | 30.52.18 | Luer bone rong str 18cm | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 32853 | GTIN | 04046826042738 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 15-243 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R001;04;ag;3V |
| Net weight (kg) | 0,189 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | MPK23-1 |



Medicon media 450.00.03 (general catalogue the art of surgery)
450.00.63 (catalogue Orthopaedics)

Complete picture

Detail images



| | | | |
|------|------|---------|--|
| size | 18cm | Coating | |
|------|------|---------|--|

Storage conditions

No special requirements

Materials

| | Material number (DIN designation) |
|-----------|-----------------------------------|
| Main part | 1.4021 (X20Cr13) |