



Product data sheet REF 30.71.03

11.12.2025

General data			
article number	30.71.03	Beyer bone rong 3mm cvd 17cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	32853	GTIN	04046826042813
Risk category	1r (EU MDR)	UMDNS-Code	15-243
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R001;04;ag;3V
Net weight (kg)	0,099	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK23-1
			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.37 (catalogue ENT) 450.00.62 (catalogue Medicon instruments Oral and Maxillofacial Surgery) 450.00.63 (catalogue Orthopaedics)		
Complete picture		Detail images	
			
size	3mm 17cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		
Screw	1.4305 (X8CrNiS18-9)		