



Product data sheet REF 31.61.36

11.12.2025

General data			
article number	31.61.36	Sypert bone rong 36cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	32853	GTIN	04046826043384
Risk category	1r (EU MDR)	UMDNS-Code	15-243
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R001;04;ag;3V
Net weight (kg)	0,350	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK23-1



Medicon media	450.00.53 (catalogue cardio)
---------------	-------------------------------

Complete picture

Detail images



size	36cm	Coating	
------	------	---------	--

Storage conditions

No special requirements

Materials

	Material number (DIN designation)
Main part Screw	1.4021 (X20Cr13) 1.4305 (X8CrNiS18-9)