



## Product data sheet REF 329.80.03

11.12.2025

| General data  |  |                                 |                         |
|---|--|---------------------------------|-------------------------|
| article number  | 329.80.03  | Chompret syndesmotom 16cm fig 3 |                         |
| Unit  | piece  | Customs code (HS code)          | 90184990                |
| GMDN  | 35559  | GTIN                            | 04046826178789          |
| Risk category   | 1r (EU MDR)  | UMDNS-Code                      | 16-667                  |
| Classification  | 6 (EU MDR)   | Base UDI-DI                     | 4046826;B;R004;04;aj;5X |
| Net weight (kg)   | 0,029  | Operating instructions          | GATDR00X                |
| Latex free  | yes  | Preparation instructions        |                         |
|   |  |                                 |                         |
| Medicon media   | 450.00.03 (general catalogue the art of surgery)<br>450.00.40 (Katalog Nr 5 Dental ) |                                 |                         |
| Complete picture  |  | Detail images                   |                         |
|    |  |                                 |                         |
| size  | 16cm   | Coating                         |                         |
| Storage conditions  |  |                                 |                         |
| No special requirements   |  |                                 |                         |
| Materials   |  |                                 |                         |
|   | Material number (DIN designation)  |                                 |                         |
| <b>Main part<br/>Handle</b>   | 1.4034 (X46Cr13)<br>1.4301 (X5CrNi18-10)   |                                 |                         |