








## Product data sheet REF 40.50.46

11.12.2025

General data			
article number	40.50.46	vacuum cleaner 230V f 40.50.40/82	
Unit	piece	Customs code (HS code)	85081100
GMDN	45590	GTIN	04046826312695
Risk category	N.A. (EU MDD)	UMDNS-Code	
Classification	NA (EU MDD)	Base UDI-DI	
Net weight (kg)	18,000	Operating instructions	
Latex free		Preparation instructions	KEIN MP
			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.63 (catalogue Orthopaedics )		
Complete picture		Detail images	
			
size	230V	Coating	
Storage conditions			
		 	
Materials			
	Material number (DIN designation)		
Main part	KEINE ANGABEN (Werkstoffangabe nicht vorhanden bzw. Stückliste)		