



Product data sheet REF 40.50.71

11.12.2025

General data			
article number	40.50.71	flat wrench pair	
Unit	box/2	Customs code (HS code)	84679900
GMDN	60909	GTIN	04046826044671
Risk category	KEIN MP (EU MDD)	UMDNS-Code	
Classification	NA (EU MDD)	Base UDI-DI	
Net weight (kg)	0,047	Operating instructions	
Latex free	yes	Preparation instructions	KEIN MP
			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.63 (catalogue Orthopaedics)		
Complete picture		Detail images	
			
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4310 (X10CrNi18-8)		