








Product data sheet REF 43.34.14

11.12.2025

| General data | | | |
|---|--|------------------------------|--------------------|
| article number | 43.34.14 | Unna comedone expressor 14cm | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 57817 | GTIN | 04046826045494 |
| Risk category | 1 (EU MDR) | UMDNS-Code | 11-084 |
| Classification | 1 (EU MDR) | Base UDI-DI | 4046826;N002;05;X3 |
| Net weight (kg) | 0,008 | Operating instructions | GATDN00X |
| Latex free | yes | Preparation instructions | |
|   | | | |
| Medicon media | 450.00.03 (general catalogue the art of surgery) | | |
| Complete picture | | Detail images | |
|    | | | |
| size | 14cm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4024 (X15Cr13) | | |