



Product data sheet REF 44.28.65

11.12.2025

| General data | | | |
|--|---|--|----------------|
| article number | 44.28.65 | sphygmomanometer aneroid w velcro cuff | |
| Unit | piece | Customs code (HS code) | 90189010 |
| GMDN | 60909 | GTIN | 04046826046378 |
| Risk category | N.A. (EU MDD) | UMDNS-Code | |
| Classification | NA (EU MDD) | Base UDI-DI | |
| Net weight (kg) | 0,505 | Operating instructions | |
| Latex free | | Preparation instructions | KEIN MP |
|  | | | |
| Medicon media | 450.00.03 (general catalogue the art of surgery) | | |
| Complete picture | | Detail images | |
|  | | | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | KEINE ANGABEN (Werkstoffangabe nicht vorhanden bzw. Stückliste) | | |
| Others | LATEX (Latex) | | |