






## Product data sheet REF 44.30.25

11.12.2025

| General data  |  |  |  |
|---|--|--|--|
| article number  | 44.30.25   | Wartenberg neurological pin wheel 18cm |  |
| Unit  | piece  | Customs code (HS code)                 | 90189084   |
| GMDN  | 47164  | GTIN                                   | 04046826046439   |
| Risk category   | 1 (EU MDR)                                       | UMDNS-Code                             | 13-043   |
| Classification  | 1 (EU MDR)                                       | Base UDI-DI                            | 4046826;N001;06;WT   |
| Net weight (kg)   | 0,057  | Operating instructions                 | GATDN00X   |
| Latex free  | yes  | Preparation instructions               |  |
|   |  |  |  |
| Medicon media   | 450.00.03 (general catalogue the art of surgery) |  |  |
| Complete picture  |  | Detail images                          |  |
|   |  |  |  |
| size  | 18cm   | Coating                                | During anodizing, titanium is enriched to form a hard metal layer on aluminum components, which can also be colored. |
| Storage conditions  |  |  |  |
| No special requirements   |  |  |  |
| Materials   |  |  |  |
|   | Material number (DIN designation)                |  |  |
| Main part   | 1.4021 (X20Cr13)                                 |  |  |
| Handle  | 3.3215 (AlSi1MgMn)                               |  |  |
| Others  | 2.0401 (CuZn39Pb3)                               |  |  |