






## Product data sheet REF 44.35.15

11.12.2025

General data			
article number	44.35.15	tape measure linen 1.5m	
Unit	piece	Customs code (HS code)	90178010
GMDN	44796	GTIN	04046826046477
Risk category	1 (EU MDR)	UMDNS-Code	
Classification	1 (EU MDR)	Base UDI-DI	4046826;N002;07;X9
Net weight (kg)	0,028	Operating instructions	GATDN00X
Latex free	yes	Preparation instructions	
 			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.63 (catalogue Orthopaedics )		
Complete picture		Detail images	
			
size	1,5m	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	ABS (Acrylnitril-Butadien-Styrol)		