






## Product data sheet REF 44.38.10

11.12.2025

General data			
article number	44.38.10	Moeltgen flexometer 23.5cm	
Unit	piece	Customs code (HS code)	90189010
GMDN	33651	GTIN	04046826046552
Risk category	1 (EU MDR)	UMDNS-Code	11-893
Classification	1 (EU MDR)	Base UDI-DI	4046826;N002;07;X9
Net weight (kg)	0,149	Operating instructions	GATDN00X
Latex free	yes	Preparation instructions	
 			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.60 (catalogue plastic surgery ) 450.00.62 (catalogue Medicon instruments Oral and Maxillofacial Surgery) 450.00.63 (catalogue Orthopaedics )		
Complete picture		Detail images	
			
size	23.5cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		