






Product data sheet REF 44.38.30

11.12.2025

General data			
article number	44.38.30	finger goniometer 14.5cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	33651	GTIN	04046826046569
Risk category	1 (EU MDR)	UMDNS-Code	11-893
Classification	1 (EU MDR)	Base UDI-DI	4046826;N002;07;X9
Net weight (kg)	0,039	Operating instructions	GATDN00X
Latex free	yes	Preparation instructions	
 			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.60 (catalogue plastic surgery) 450.00.63 (catalogue Orthopaedics)		
Complete picture		Detail images	
			
size	14.5cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4301 (X5CrNi18-10)		