




Product data sheet REF 44.51.70

11.12.2025

General data			
article number	44.51.70	diagnostic station 230V oto/ophthalmoscope	
Unit	piece	Customs code (HS code)	90189084
GMDN	60909	GTIN	04046826046897
Risk category	N.A. (EU MDD)	UMDNS-Code	
Classification	NA (EU MDD)	Base UDI-DI	
Net weight (kg)	3,564	Operating instructions	
Latex free		Preparation instructions	KEIN MP
			
Medicon media			
Complete picture		Detail images	
size	230V	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	KEINE ANGABEN (Werkstoffangabe nicht vorhanden bzw. Stückliste)		