





## Product data sheet REF 45.30.08

11.12.2025

| General data  |  |                                |                    |
|---|--|--------------------------------|--------------------|
| article number  | 45.30.08                                       | Miller f.o. laryng blade fig 3 |                    |
| Unit  | piece  | Customs code (HS code)         | 90189084           |
| GMDN  | 46827  | GTIN                           | 04046826535827     |
| Risk category   | 1 (EU MDR)                                     | UMDNS-Code                     | 15-076             |
| Classification  | 10 (EU MDR)                                    | Base UDI-DI                    | 4046826;N004;05;XR |
| Net weight (kg)   | 0,000  | Operating instructions         | GATDN00X           |
| Latex free  | yes  | Preparation instructions       |                    |
|   |  |                                |                    |
| Medicon media   | 453.11.01 (Fiber Optic Laryngoscopes Brochure) |                                |                    |
| Complete picture  |  | Detail images                  |                    |
|   |  |                                |                    |
| Storage conditions  |  |                                |                    |
| No special requirements   |  |                                |                    |
| Materials   |  |                                |                    |
|   | Material number (DIN designation)              |                                |                    |
| Main part   | 1.4301 (X5CrNi18-10)                           |                                |                    |