





## Product data sheet REF 45.30.15

11.12.2025

General data			
article number	45.30.15	xenon lamp 2.5V only	
Unit	piece	Customs code (HS code)	90189084
GMDN	60909	GTIN	04046826535407
Risk category	N.A. (EU MDR)	UMDNS-Code	
Classification	NA (EU MDR)	Base UDI-DI	
Net weight (kg)	0,000	Operating instructions	GATDN004
Latex free	yes	Preparation instructions	
 			
Medicon media	453.11.01 (Fiber Optic Laryngoscopes Brochure)		
Complete picture		Detail images	
size	2.5V	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		