



Product data sheet REF 45.30.16

11.12.2025

General data			
article number	45.30.16	xenon lamp 3.5V only	
Unit	piece	Customs code (HS code)	90189084
GMDN	60909	GTIN	04046826535070
Risk category	N.A. (EU MDR)	UMDNS-Code	
Classification	NA (EU MDR)	Base UDI-DI	
Net weight (kg)	0,000	Operating instructions	GATDN004
Latex free	yes	Preparation instructions	



Medicon media 453.11.01 (Fiber Optic Laryngoscopes Brochure)
453.32.00 (FlexiTrac - Abdominal retractor system Brochure)

Complete picture

Detail images

size	3.5V	Coating	
------	------	---------	--

Storage conditions

No special requirements

Materials

Material number (DIN designation)