



Product data sheet REF 45.30.18

11.12.2025

General data			
article number	45.30.18	hdl slim w inductive accu w xenon lamp	
Unit	piece	Customs code (HS code)	90189084
GMDN	46830	GTIN	04046826534400
Risk category	1 (EU MDR)	UMDNS-Code	15-076
Classification	10 (EU MDR)	Base UDI-DI	4046826;N001;06;WT
Net weight (kg)	0,000	Operating instructions	GATDN00X
Latex free	yes	Preparation instructions	



Medicon media
 453.11.01 (Fiber Optic Laryngoscopes Brochure)
 453.32.00 (FlexiTrac - Abdominal retractor system Brochure)

Complete picture

Detail images



Storage conditions

No special requirements

Materials

	Material number (DIN designation)
Main part	1.4301 (X5CrNi18-10)